

## - COL G -

### Recombinant Collagenase class I



Item No.	Item Description
001-0025	COL G, 25 mg ≥ 75 U
001-0100	COL G, 100 mg ≥ 300 U
001-0250	COL G, 250 mg ≥ 750 U

## - COL H -

### Recombinant Collagenase class II



Item No.	Item Description
002-0025	COL H, 25 mg ≥ 750 U
002-0100	COL H, 100 mg ≥ 3000 U
002-0250	COL H, 250 mg ≥ 7500 U

#### 1. DESCRIPTION

**COL G** and **COL H** are recombinant collagenases (metalloproteinases) class I and class II respectively [1]. COL G and COL H are synthesized separately from *C. Histolyticum* genes by DNA recombination in *E. Coli BL21 AI* strain, bearing a Maltose Binding Protein (MBP) tag at the N-terminal end [2].

COL G and COL H are affinity chromatography purified proteins, highly pure, highly stable, lot-to-lot consistent, endotoxin-free ( $\leq 10$  EU/mg, LAL assay) and animal-free.

<b>CAS:</b>	9001-12-1
<b>EC:</b>	3.4.24.3
<b>Grade:</b>	Research Premium Grade
<b>Form:</b>	Lyophilized white powder
<b>Quality:</b>	Amylose Affinity Chromatography
<b>Inhibitors:</b>	EDTA, EGTA, Cys, Hys, DTT, 2-mercaptoethanol
<b>Activators:</b>	Ca <sup>2+</sup>

Their molecular weights are ~135 kDa (COL G) and ~158.5 kDa (COL H). COL G and COL H are soluble in water or aqueous buffers and express their maximum activity at **pH 8**.

#### 2. SUBSTRATES

COL G and COL H play different synergic roles in collagen digestion. Indeed, COL G expresses a higher activity against **native collagen**, specifically hydrolyzing **3D-helix regions**, while

COL H expresses a lower activity against the 3D helix and a higher activity against **linear collagen regions** at the motif Pro-Y-Gly-Pro [3,4]. The mix of COL G and COL H expresses a **synergic activity** that results in efficient collagen digestion [5].

For tissue dissociation, protease addition is needed to hydrolyze non-collagenous proteins and other macromolecules present in the extracellular matrix [6].

#### 3. ENZYMATIC ACTIVITY

**COL G** ≥ 3.0 Units/mg\*  
**COL H** ≥ 30.0 Units/mg\*

\*according to Grassmann, one Unit liberates 1  $\mu$ mol of Gly-Pro-Ala from Carbobenzoxy-Gly-Pro-Gly-Gly-Pro-Ala-OH (Fluka 27673) in 1 min at pH 7.4, 37 °C [7].

#### 4. APPLICATIONS

**For research use only.**

Due to their high purity and specificity, COL G and COL H are especially indicated for the isolation of primary cells from liver, pancreas, heart, cartilage and stem cells from adipose tissue and others.

In these applications we recommend using a combination of COL G and COL H in a specific activity ratio, or according to the relevant isolation protocol in order to obtain an optimal collagen digestion in cell isolation. For other applications or suggestions, contact [info@abielbiotech.com](mailto:info@abielbiotech.com) or visit [www.abielbiotech.com](http://www.abielbiotech.com).

### 5. PREPARATION METHOD

We recommend reconstituting the lyophilized COL G and COL H enzymes in the tissue-dissociation buffer by injecting the **buffer directly into the vial**. Do not exceed an enzyme concentration of 30 U/ml (COL G) or 300 U/ml (COL H) to avoid precipitates.

Keep the vial on ice and periodically shake until the enzyme is completely dissolved. Filter with 0.22 µm mesh for sterility.

**Prepare a mix of COL G and COL H** solutions in a specific activity ratio and dilute according to your protocol working solution concentration.

Add protease to the mix at 4 °C according to the specific application. Thermolysin, pronase or neutral protease/dispase can be normally used.

**Protease must be added immediately before use** to avoid catalytic processes in the enzymatic blend. The amount of protease will define the aggressiveness of your enzyme mixture. For suggestions about your specific protocol and application please contact [info@abelbiotech.com](mailto:info@abelbiotech.com) or visit [www.abelbiotech.com](http://www.abelbiotech.com).

### 6. STORAGE AND STABILITY

Lyophilized COL G and COL H are stable at -80 °C up to two years. We recommend splitting in aliquots the reconstituted solutions at need and storing them at -20°C up to one month or -80°C up to 6 months.

To use aliquots later on, they can be diluted in re-constitutive buffer or can be directly added into the enzyme working solution.

**▲ Warning: We recommend avoiding multiple freeze-thaw cycles and exposure to frequent temperature changes.**

#### REFERENCES

- [1] Matsushita, O. et al. (1999) *J. Bacteriol.* 181(3): 923–933.
- [2] Salamone, M. et al. (2012) *Chem. Eng. Trans.* 27: 259-264.
- [3] Philominathan S.T. et al. (2009) *J. Biol. Chem.* 284(16) : 10868-10876;
- [4] Matsushita, O. Et al. (1994) *J. Bacteriol.* 176: 149-156
- [5] Breite, A.G. et al. (2011) *Transplant Proc.* 43(9) : 3171-3175
- [6] Salamone, M. et al. (2014) *Chem. Eng. Trans.* 38: 247-252.
- [7] W. Grassmann, et al, (1960) *Z. Physiol.Chemie* 322:267

**For suggestions about your specific protocol or application of COL G and COL H, contact us:**

Monica Salamone, PhD  
R&D Manager  
[m.salamone@abelbiotech.com](mailto:m.salamone@abelbiotech.com)

Alessandra Lo Cicero, PhD  
Product Manager  
[a.locicero@abelbiotech.com](mailto:a.locicero@abelbiotech.com)

Tel. +39 091 6615633  
Fax +39 091 661 5628